


ARTERY CLAMPS AND GOUGE FORCEPS



American Eagle clamps and gouge pliers - manufactured to the highest standards of surgical mechanical craftsmanship and maximum precision.

Fine, precise cutting edges enable the gouge pliers to be **used precisely and smoothly**. A standard version as well as a particularly fine variant are available depending on the clinical indication. For example, bones can be contoured, but fine bone chips can also be carved out of a block.

The finely crafted clamps have a particularly precise instrument closure. This ensures that the fine knurlings on the working ends close properly to stop bleeding and minimize tissue trauma.



AESCHMS	Halsted-Mosquito Haemostatic Forceps, Straight, 14 cm Halsted-Mosquito haemostatic forceps are suitable for clamping small blood vessels.
AESCHMC	Halsted-Mosquito Haemostatic Forceps, Angled, 14 cm Halsted-Mosquito haemostatic forceps are suitable for clamping small blood vessels.
AESCMMOS	Micro Mosquito Haemostatic Forceps, Straight, 12 cm Micro Mosquito haemostatic forceps are suitable for clamping small blood vessels.
AESCMMOC	Micro Mosquito Haemostatic Forceps, Angled, 12 cm Micro Mosquito haemostatic forceps are suitable for clamping small blood vessels.

AESBRFRIE	Friedmann Bone Rongeur Forceps Friedmann bone rongeur forceps are used to contour bone and thick tissue.
AESBRFRIEM	Mini-Friedmann Bone Rongeur Forceps The particularly fine Mini-Friedmann bone rongeur forceps are used to contour bone and thick tissue.