




DIAGNOSTICS AND ACCESSORIES





In addition to the surgical instruments from American Eagle Instruments®, diagnostic instruments and accessory are also available in a matching design.

To avoid disturbing light reflections, the handles are matted. The integrated macro-grooves provide a **good grip** and are **easy to clean** at the same time.

AESH51	Dental Mirror Handle The ergonomically shaped mirror handle is suitable for commercially available dental mirrors with a metric thread (European thread).
AESH51E	Dental mirror handle with endodontic measuring scale (0-30 mm) The ergonomically shaped mirror handle is suitable for commercially available dental mirrors with a metric thread (European thread). With an easy-to-read endodontic measuring scale from 0 to 30 mm on the opposite end of the handle.
AESPUNC15	UNC 15 Periodontal Probe, University of North Carolina, 1-15 mm, 1 mm Scale Increments UNC-15 Periodontal probe with 1 mm scale increments. Suitable for simple measurements/assessments, such as implant spacing, crown lengths, bone thicknesses, etc.
AESCT10	London College Tweezers, smooth handle London College multi-purpose tweezers with smooth handle. Particularly easy to clean and sterilise. Suitable for a wide variety of procedures.
AESCT5	London College Tweezers, serrated handle London College Multi-Purpose Tweezers with Serrated Handle. With guide pin for maximum hold and a serrated handle for a secure grip. Suitable for a wide variety of procedures.
AESGP	Luniatschek Gauze Packer For the easy placement of gauze strips/tamponades. The matt handle surface reduces irritating light reflections.