

GALAXIE
cassettes

GALAXIE INSTRUMENT CASSETTES


The cost efficient solution for instrument reprocessing

NEW!


Innovative tooth design
Additional cassette size
New accessories



GALAXIE cassettes



innovative
tooth design



AEC1/4DIN

Cassette for 5 instruments, ¼ DIN, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black), 185 x 80 x 37 mm

Extendable to a 10-instrument cassette with additional Retention Rail (1x AECAGRAIL)

The gold-standard solution for instrument reprocessing

The new instrument cassettes with innovative tooth design guarantee maximum efficiency in the cleaning and sterilisation process. Additional cassette models, such as container and 1/1 DIN surgery cassette, as well as new clips for attaching hinged instruments increase the range of options for optimum reprocessing organisation in your practice.

High-quality Materials:

- » Electropolished stainless steel makes the cassette edges smooth and protects against corrosion
- » Only food grade silicone used

Innovative and ergonomic closing latch:

- » Safe and easy-to-use locking system, allows one-handed opening


New design for a better cleaning process:

- » The new and innovative laser cutting pattern guarantees maximum efficiency during the cleaning and sterilization process whilst ensuring stability and strength of the cassette
- » Open lattice design specially designed for washer disinfectant use, enabling maximum water flow
- » Smooth round corners and slotted edges increase drainage and reduce drying time


GALAXIE cassettes

AEC1/2DIN

Cassette for 10 instruments, ½ DIN, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black), 185 x 148 x 37 mm




Colour coding
option



Innovative color coding and modular silicone holders:

- » The innovative shape reduces instrument contact allowing for more liquid flow and for more efficient cleaning
- » 5 color coding tabs for practice organization included with every cassette
- » Silicone rails can be assembled in any position, are compatible with all handle diameters

Extendable to a 20-instrument cassette with additional Retention Rails (2x AECAGRAIL)



GALAXIE cassettes


AESCCON

Container cassette for up to 16 instruments, fits DIN dental standard containers, with handle on lid, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black), 268 x 176 x 34 mm



GALAXIE cassettes

Accessories



AECRAIL

Standard instrument rail, blue for 10 instruments, 1 pc. Suitable for the GALAXIE INSTRUMENT CASSETTES AEC1/2DIN, AECDIN, AESCDIN and AESCCON.

For use in the AEC1/4DIN cassette, the rail can be cut in half.

AECAGRAIL

Retention instrument rail, blue for 10 instruments, 1 pc. Suitable for the GALAXIE INSTRUMENT CASSETTES AEC1/2DIN, AECDIN, AESCDIN and AESCCON.

For use in the AEC1/4DIN cassette, the rail can be cut in half.



AECRAIL

Compression rail/ pressure rail V-profile, blue for 10 instruments, 1 pc. Suitable for the GALAXIE INSTRUMENT CASSETTES AEC1/2DIN, AECDIN, AESCDIN and AESCCON.

For use in the AEC1/4DIN cassette, the rail can be cut in half.

AESCLIPNH

Clamp for needle holders or forceps, 1 pc. Secures needle holders/clamps/forceps in the lid of the instrument cassette. NOTE: Use in pairs recommended for optimal fixation.

Suitable for all GALAXIE Instrument Cassettes.



AESCLIPS

Clamp for scissors, 1 pc. Secures scissors in the lid of the instrument cassette. NOTE: Use in pairs recommended for optimal fixation.

Suitable for all GALAXIE Instrument Cassettes.



AESCTOOL

Assembly tool for scissors and needle holder clips, 1 pc.

Suitable for all GALAXIE Instrument Cassettes.




GALAXIE cassettes

Item #	Description
AEC1/4DIN	Cassette for 5 instruments, ¼ DIN, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black) 185 x 80 x 37 mm, with additional Retention Rail (1x AECAGRAIL) expandable to 10 instruments
AEC1/2DIN	Cassette for 10 instruments, ½ DIN, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black) 185 x 148 x 37 mm, with additional Retention Rails (2x AECAGRAIL) expandable to 20 instruments
AESCDIN	Surgical cassette for up to 17 instruments, DIN, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black), 285 x 183 x 37 mm
AESCCON	Container cassette for up to 16 instruments, with handle on lid, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black), 268 x 176 x 34 mm, fits standard DIN sterilisation containers
AECAGRAIL	10-instruments Retention Rail, blue, 1 piece
AECRAIL	10-instruments Support Rail, blue, 1 piece
AECRAIL	Pressure Rail, blue, V-profile, 1 piece
AESCLIPS	Clip for scissors, 1 piece (NOTE: For optimal fixation, use in pairs is recommended)
AESCLIPNH	Clip for needle holders or pliers, 1 piece (NOTE: Use in pairs is recommended for optimal fixation)
AESCTOOL	Assembly tool for scissors and needle holder clips, 1 piece

Cassettes are supplied without instruments and without additional silicone rails.

AESCDIN

Surgical cassette for up to 17 instruments, DIN, blue rails, incl. 5 colour coding options (purple, red, yellow, green, black), 285 x 183 x 37 mm



AEI-GB-CAS/1123