



PERIOLUX ELEVATORS



Developed for the gentlest possible extraction. The instruments were designed for axial use in the pencil grip which enables direct and controlled power transmission for maximum control. Surrounding bone structures can be protected as much as possible and tooth / root fractures can be avoided.

The matt handles reduce irritating light reflection and are ergonomically shaped, which also ensures direct and controlled power transmission. The MacroGrip grooved structure of the handles offers a secure hold and is easy to clean.

PerioLux instruments are not suitable for levering!

AESPLXD	PerioLUX Distal, ergonomic, for distal use
AESPLXM	PerioLUX Mesial, ergonomic, for mesial use

