




# PERIOTOMES



The American Eagle periotomes were developed for the gentlest possible extraction, especially before planned implantation. Before using dental forceps or root elevators, the periodontal ligament can be separated with periotomes as **atraumatic as possible**. The very fine working ends are designed for this purpose, to be able to penetrate deep into the periodontal gap.

The hardened, fine and slightly flexible instrument blades follow the root structure and separate the fibres of the periodontium. The instruments were designed with a pencil grip for maximum control. They enable **direct and controlled power transmission** and the surrounding bone structures are protected as much as possible. The matte handles avoid annoying light reflections.

The MacroGrip grooved structure of the handles offers a secure hold and is easy to clean.

<b>AESPT4</b>	Periotome 1,8mm, staight and curved
<b>AESPT5</b>	Periotome 2,5mm straight and curved
<b>AESPT6</b>	Periotome 1,8mm and 2,5mm, curved

