


RETRACTOR / CHEEK RETRACTOR





American Eagle retractors are made of surgical steel for optimal cleaning and preparation. The **mat surface finish** avoids annoying light reflections. The models with a retaining ring have been specially modified for processing in the appropriate American Eagle surgical cassettes - due to the **adjusted ring size**. These can be processed together with the other instruments in the closed cassette.

AESRCW	Cawood-Minnesota Retractor The Cawood-Minnesota retractor was specially designed to meet a variety of needs in dental surgery: for retracting mucoperiosteal flaps, the cheeks, lips or tongue.
AESRLAN1	Langenbeck Cheek Retractor 1, 30 x 16 x 210 mm
AESRLAN2	Langenbeck Cheek Retractor 2, 40 x 11 x 210 mm
AESRLAN3	Langenbeck Cheek Retractor 3, 50 x 11 x 210 mm The Langenbeck cheek retractor 3 is used to retract mucoperiosteal flaps, the cheeks, lips or tongue.
AESRMID1	Middeldorpf Retractor #1, 14 x 17 x 210 mm
AESRMID2	Middeldorpf Retractor #2, 20 x 22 x 210 mm The Middeldorpf cheek retractor 2 is used to retract mucoperiosteal flaps, the cheeks, lips or tongue.
AESRVES	Branemark Vestibulum Retractor The Branemark retractor, also known as an implant retractor, ensures that the surgical area is easily accessible and clear.
AESRB	Bishop Cheek Retractor The Bishop retractor was specially designed to meet a wide range of requirements in dental surgery and optimally follows the oral anatomy- for retracting mucoperiosteal flaps, the cheeks, lips or tongue.