



ROOT ELEVATORS



American Eagle LiftaLUX elevators - developed for **gentle extraction**.

The careful use of a delicate elevator helps to master and overcome the forces that fix the tooth in the best possible way.

The fine instrument tips should be applied axially with controlled pressure, preferably around the entire tooth, to evenly expand the socket. The ergonomic handles lie well in the hand and enable **controlled power transmission**.

LiftaLUX instruments are not suitable for levering!

AESRES3	LiftaLUX Elevator, straight, ergonomic, 3mm 3 mm straight working end. LiftaLUX instruments are not suitable for levering!
AESRES4	LiftaLUX Elevator, straight, ergonomic, 4mm 4 mm straight working end. LiftaLUX instruments are not suitable for levering!
AESREC3	LiftaLUX Elevator, curved, ergonomic, 3mm 3 mm curved working end for better lateral access. LiftaLUX instruments are not suitable for levering!
AESREC4	LiftaLUX Elevator, curved, ergonomic, 4mm 4 mm curved working end for improved lateral access. LiftaLUX instruments are not suitable for levering!



TwistaLUX instruments are suitable for increased force and **levering and torque extraction**. The working ends are more robust here.



AESRETS3	TwistaLUX Elevator, straight, ergonomic, 3mm 3 mm straight working end - specially designed for levering and torque applications.
AESRETS4	TwistaLUX Elevator, straight, ergonomic, 4mm 4 mm straight working end - specially designed for levering and torque applications.
AESRETC3	TwistaLUX Elevator, curved, ergonomic, 3mm 3 mm curved working end for better lateral access - specially designed for levering and torque applications.
AESRETC4	TwistaLUX Elevator, curved, ergonomic, 4mm 4 mm curved working end for better lateral access - specially designed for levering and torque applications.