




ROOT TIP ELEVATORS




YOUNG[®]
INNOVATIONS



The Root tip lifter removes the smallest root fragments. They are manufactured with an ergonomic handle for **controlled power transmission**.

The instrument tips are hardened and lengthened for better/deeper access to the socket. The **matt handle surface** reduces annoying light reflections.



AESREHEID1	Heidbrink Root Tip Elevator #1, fine, ergonomic
AESREHEID2	Heidbrink Root TipElevator #2, medium, ergonomic
AESREHEID3	Heidbrink Root Tip Elevator #3, strong, ergonomic